



ADDITIONAL INFORMATION

Weight	0.04 kg
Dimensions	9.7 × 3.5 × 11.6 cm
Carton Qty	20
Pricing Unit	Box/10
Size	5 x 5cm
Min Order Qty	1
ARTG ID	406570
GMDN Code	43186
Medical Device Class	Class 2b
Tax Schedule	Non-Taxable
Item Length	5
Item Width	5
Item Measurement Unit	cm
Is Medical Device	Yes
Is Serialized Item	No
Is Dangerous Good	No



PRODUCT DESCRIPTION

AeroWound Alginate Dressings provide an advanced wound care solution, designed to provide a moist healing environment. Naturally derived from seaweed, AeroWound is an ideal choice for various wound types. On contact with wound exudate, the dressing transforms into a soft gel, maintaining a moist healing environment which is conducive to wound healing, and assists in pain reduction.

<https://aerohealthcareonline.com/wp-content/uploads/2023/04/AlginateDressing-Final.mp4>

Features:

1. **Optimized Healing Environment:** Our dressings ensure a moist wound environment, speeding up the natural healing process.
2. **High Absorbency:** Specially crafted to absorb wound exudate efficiently, minimising the risk of maceration.
3. **Natural Seaweed Derivative:** Made from natural seaweed, providing bio-degradable, eco-friendly care.
4. **Soft & Conforming:** Flexible design adheres comfortably to the wound contours, ensuring patient comfort.
5. **Easy Removal:** Designed to gel up upon contact with exudate for painless removal, reducing patient trauma.

Safety Considerations: For external use only. Always read the label and follow the directions for use.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Alginate Dressing 5 x 5cm