AEROWOUND COMBINE ROLL 10CM X 10M



ACD1010R



ADDITIONAL INFORMATION

Weight 0.268 kg

Dimensions $12 \times 12 \times 10 \text{ cm}$

Carton Qty 10

Pricing Unit Dressing

Size 10cm x 10M

Min Order Qty 1

ARTG ID 404994
GMDN Code 45050
Medical Device Class 1
Tax Schedule Taxable

Item Length 1000

Item Width 10

Item Measurement Unit cm

Is Medical Device Yes

Is Serialized Item No

Is Dangerous Good No

PRODUCT DESCRIPTION

The AeroWound Combine Dressing is recommended where a highly absorbent dressing is required for large and heavily bleeding wounds.

Low Adherent


• The AeroWound Combine Dressing has a special low adherent coating which minimizes trauma associated with changing wound dressings

Hypo-Allergenic



AEROWOUND COMBINE ROLL 10CM X 10M

ACD1010R



Highly Absorbent

• The AeroWound Combine Dressing is capable of absorbing fluids from large and heavily bleeding wounds to stop bleeding and reduce infection risk.

Individually packed and sterile

• Dressing packs are available individually on demand in an easy peel sterile pouch to ensure wound cleanliness.

Range of Sizes Available

• Ranges of dressing sizes are available to suit both small and large sized wounds effectively.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Combine Roll 10cm x 10M