

AEROWOUND COMBINE DRESSING 10 X 10CM BAG/50

ACD1010S





ADDITIONAL INFORMATION

Weight 0.27 kg

Dimensions $32 \times 26 \times 2 \text{ cm}$

Carton Qty 16

Pricing Unit Bag/50

Size 10cm x 10cm

Min Order Qty 1

ARTG ID 159437 **GMDN Code** 34654

Medical Device Class 1s

UNSPSC 42311545

Tax Schedule Taxable

Item Length 10
Item Width 10

Item Measurement Unit cm

Is Serialized Item No

Is Dangerous Good No

PRODUCT DESCRIPTION

The AeroWound Combine Dressing is recommended where a highly absorbent dressing is required for large and heavily bleeding wounds.

Low Adherent

• The AeroWound Combine Dressing has a special low adherent coating which minimizes trauma associated with changing wound dressings



AEROWOUND COMBINE DRESSING 10 X 10CM BAG/50

ACD1010S

Hypo-Allergenic

 The Aer 	oWound Com	bine Dressing	is hypo-allerg	enic with a	low risk of we	ound irritation
-----------------------------	------------	---------------	----------------	-------------	----------------	-----------------

Highly Absorbent

• The AeroWound Combine Dressing is capable of absorbing fluids from large and heavily bleeding wounds to stop bleeding and reduce infection risk.

Individually packed and sterile

• Dressing packs are available individually on demand in an easy peel sterile pouch to ensure wound cleanliness.

Range of Sizes Available

• Ranges of dressing sizes are available to suit both small and large sized wounds effectively.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Combine Dressing bag/50