

## AEROWOUND COMBINE DRESSING 20 X 40CM

**ACD2040S** 





## **ADDITIONAL INFORMATION**

Weight 0.03 kg

**Dimensions** 28 × 15 × 1 cm

**Carton Qty** 100

**Pricing Unit** Unit/1

Size 20cm x 40cm

Min Order Qty

**ARTG ID** 159437 **GMDN Code** 34654

**Medical Device Class** Class 1s

UNSPSC 42311545

**Tax Schedule** Taxable

Item Length 40

**Item Width** 20

Item Measurement Unit cm

Is Medical Device Yes

Is Serialized Item Nο

Is Dangerous Good No



## PRODUCT DESCRIPTION

AeroWound Combine Dressings provide high-capacity wound padding for situations where swift absorption and dependable cushioning are essential. Built for post-surgical sites, burns and heavily bleeding wounds, the soft, cotton-filled core draws and holds exudate while helping protect the area from further trauma. A low-adherent contact layer is designed to minimise sticking on removal, supporting comfortable dressing changes and continuity of care across clinical and first aid environments.

The dressing utilises a nonwoven outer wrap with generous loft to deliver bulk absorption without excessive weight. Each piece is sterile (ethylene oxide) and intended for single use only. Individual, easy-peel pouches support aseptic presentation at the point of care and help maintain pack integrity in kits, trolleys and emergency response bags. Clear on-pack icons communicate key handling and storage requirements, including temperature limits and 'do not re-use' guidance.

Combine Dressings integrate seamlessly into wound-care workflows. Apply directly over the wound as a primary pad or use as a secondary absorbent layer over topical treatments or contact dressings. The square format sits neatly under retention bandages and fixation tapes, making it a practical option for hospitals, day clinics, aged care, workplace first aid stations and community responders.



## **AEROWOUND COMBINE DRESSING 20 X 40CM**

**ACD2040S** 

- **High absorbency:** cotton-filled core manages heavy exudate
- Low-adherent contact layer: supports gentler removal
- Sterile, single use: ethylene oxide sterilised, pouch-packed
- **Nonwoven wrap:** soft, low-lint cushioning for protection
- Easy handling: individual easy-peel sterile pouch
- Clear storage guidance: keep dry; 0-30°C temperature limit
- Versatile: suitable as primary or secondary absorbent pad

The **10 × 10 cm** option delivers compact, square coverage with ample absorbency for post-surgical sites and burns, fitting neatly under fixation tapes and bandages for secure retention.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Combine Dressing dressing/1