



ADDITIONAL INFORMATION

Weight	0.18 kg
Dimensions	27.1 × 4.3 × 29.6 cm
Carton Qty	20
Pricing Unit	Bag/50
Size	4.5cm x 4.5cm
Min Order Qty	1
ARTG ID	159433
GMDN Code	34970
Medical Device Class	Class 1s
UNSPSC Code	42311545
Tax Schedule	Taxable
Is Medical Device	Yes
Is Serialized Item	No
Is Dangerous Good	No
Shelf Life (Months)	60
Warranty Period (Months)	0



PRODUCT DESCRIPTION

AEROWOUND Finger Dressing is an all-in-one solution designed specifically for finger injuries that are difficult to cover with standard adhesive strips. It combines a 4.5 cm low-adherent, absorbent pad with an integrated gauze wrapping leader and a secure adhesive tab to create focused pressure and reliable coverage. The construction helps control bleeding and shields the wound from everyday contaminants. Each unit is sterile and individually enclosed for cleanliness and confidence at the point of use.

The soft, low-adherent contact layer minimises trauma during removal while the absorbent pad manages light exudate for comfortable wear between changes. A superior-stretch wrap conforms closely to finger contours and joints, maintaining a neat, secure

hold during movement at work, home, or on the go. The dressing is latex-free and single use, making it a practical choice for first aid kits, workplace safety stations, mobile responders and training environments where a quick, hygienic response is essential.

By combining pad, wrap and closure in one sterile unit, AEROWOUND Finger Dressing reduces the number of separate components needed in kits and speeds up application when seconds matter. The integrated adhesive tab eliminates the need for extra fixation tape for most applications, simplifying set-up and helping users achieve a tidy, secure finish.

- **All-in-one convenience:** pad, wrap and adhesive tab simplify finger dressing.
- **Sterile and single use:** hygienic application when it matters.
- **Low-adherent pad:** helps reduce trauma on removal.
- **Superior-stretch wrap:** conforms to fingers for a secure fit.
- **Compact pouch:** easy to stock in first aid kits and vehicles.

TEst

Supplied in a **bag of 50** sterile, single-use finger dressings for efficient kit replenishment in workplaces, education and field teams. Bulk format supports fast response and fewer restocks while keeping storage compact.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Finger Dressing