





## **ADDITIONAL INFORMATION**

Weight 0.18 kg

**Dimensions**  $27.1 \times 4.3 \times 29.6 \text{ cm}$ 

Carton Qty 20

Pricing Unit Bag/50

**Size** 4.5cm x 4.5cm

Min Order Qty 1

**ARTG ID** 159433 **GMDN Code** 34970

Medical Device Class 1s

**UNSPSC** 42311545

Tax Schedule Taxable

Is Serialized Item No
Is Dangerous Good No
Shelf Life (Months) 60

### **PRODUCT DESCRIPTION**

The AEROWOUND Wound Dressings are a general purpose wound dressing for use in preventing infection and controlling bleeding. Available in sizes to suit wounds from finger injuries to large wounds.

# **Highly Absorbent**

Highly absorbent wound pad controls bleeding and absorbs wound exudates

## Low adherent

Special low adherent coating minimizes wound trauma associated with changing dressings



# **AEROWOUND FINGER DRESSING BAG/50**

AWF002

## Individually packed and sterile

• Dressings are available individually on demand in an easy tear wrapper and sterile for wound cleanliness.

# **Ease of application**

• The wound pad is attached to a wrapping leader of gauze bandage making the dressing simple to apply, yet comfortable and secure for all general wound dressing.

### **Unique Finger Dressing**

• The AEROWOUND Finger dressing is unique in the industry for it ability to treat finger wounds too large to cover with an adhesive dressing. This unique product combines a 4.5cm square low adherent and absorbent pad, with a wrapping leader and secure adhesive tab to protect the wound.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Finger Dressing