





ADDITIONAL INFORMATION

Weight 0.18 kg

Dimensions $27.1 \times 4.3 \times 29.6 \text{ cm}$

Carton Qty 20

Pricing Unit Bag/50

Size 4.5cm x 4.5cm

Min Order Qty 1

ARTG ID 159433 **GMDN Code** 34970

Medical Device Class 1s

UNSPSC 42311545

Tax Schedule Taxable

Is Medical Device Yes
Is Serialized Item No
Is Dangerous Good No
Shelf Life (Months) 60

PRODUCT DESCRIPTION

The AEROWOUND Wound Dressings are a general purpose wound dressing for use in preventing infection and controlling bleeding. Available in sizes to suit wounds from finger injuries to large wounds.

Highly Absorbent

• Highly absorbent wound pad controls bleeding and absorbs wound exudates

Low adherent

· Special low adherent coating minimizes wound trauma associated with changing dressings

Individually packed and sterile

• Dressings are available individually on demand in an easy tear wrapper and sterile for wound cleanliness.



AEROWOUND FINGER DRESSING BAG/50

AWF002

Ease of application

• The wound pad is attached to a wrapping leader of gauze bandage making the dressing simple to apply, yet comfortable and secure for all general wound dressing.

Unique Finger Dressing

• The AEROWOUND Finger dressing is unique in the industry for it ability to treat finger wounds too large to cover with an adhesive dressing. This unique product combines a 4.5cm square low adherent and absorbent pad, with a wrapping leader and secure adhesive tab to protect the wound.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Finger Dressing