AEROWOUND FINGER DRESSING ENV/3







ADDITIONAL INFORMATION

Weight 0.02 kg

Dimensions $15.3 \times 8.5 \times 2.6 \text{ cm}$

Carton Qty 96

Pricing Unit Env/3

Size 4.5cm x 4.5cm

Min Order Qty 1

ARTG ID 159433 **GMDN Code** 34970

Medical Device Class 1s

UNSPSC 42311545

Tax Schedule Taxable

Item Length 4.5

Item Width 4.5

Item Measurement Unit cm

Is Medical Device Yes

Is Serialized Item No

Is Dangerous Good No

Shelf Life (Months) 60

PRODUCT DESCRIPTION

AEROWOUND Finger Dressing is an all-in-one solution designed specifically for finger injuries that are difficult to cover with standard adhesive strips. It combines a 4.5 cm low-adherent, absorbent pad with an integrated gauze wrapping leader and a secure adhesive tab. The design enables focused pressure and reliable coverage to help control bleeding and protect the wound from everyday contaminants. Each unit is sterile and individually enclosed for cleanliness and confidence at the point of use.

The soft, low-adherent contact layer minimises trauma during removal while the absorbent pad manages light exudate, supporting comfortable wear between changes. A superior-stretch wrap conforms closely to finger contours and joints, maintaining a neat, secure hold during movement. The dressing is latex-free and single use, making it suitable for workplaces, training environments, first aid kits, and mobile responders who require a practical, hygienic option for finger-specific care.

• All-in-one convenience: pad, wrap, and adhesive tab simplify finger dressing.



AEROWOUND FINGER DRESSING ENV/3

AWF03

- Sterile and single use: hygienic application when it matters.
- Low-adherent pad: helps reduce trauma on removal.
- **Superior stretch wrap:** conforms to fingers for a secure fit.
- **Compact unit pack:** easy to stock in first aid kits and vehicles.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Finger Dressing