



ADDITIONAL INFORMATION

Weight	0.02 kg
Dimensions	15.3 × 8.5 × 2.6 cm
Carton Qty	96
Pricing Unit	Env/3
Size	4.5cm x 4.5cm
Min Order Qty	1
ARTG ID	159433
GMDN Code	34970
Medical Device Class	Class 1s
UNSPSC	42311545
Tax Schedule	Taxable
Item Length	4.5
Item Width	4.5
Item Measurement Unit	cm
Is Serialized Item	No
Is Dangerous Good	No
Barcode (GTIN/EAN/UPC)	9341394005259

PRODUCT DESCRIPTION

The AEROWOUND Wound Dressings are a general purpose wound dressing for use in preventing infection and controlling bleeding. Available in sizes to suit wounds from finger injuries to large wounds.

Unique Finger Dressing

- The AEROWOUND Finger dressing is unique in the industry for its ability to treat finger wounds too large to cover with an adhesive dressing. This unique product combines a 4.5cm square low adherent and absorbent pad, with a wrapping leader and secure adhesive tab to protect the wound.

Highly Absorbent

- Highly absorbent wound pad controls bleeding and absorbs wound exudates

Low adherent

- Special low adherent coating minimizes wound trauma associated with changing dressings

Individually packed and sterile

- Dressings are available individually on demand in an easy tear wrapper and sterile for wound cleanliness.

Ease of application

- The wound pad is attached to a wrapping leader of gauze bandage making the dressing simple to apply, yet comfortable and secure for all general wound dressing.

Finger Dressing