

AEROWOUND #13 WOUND DRESSING 8 X 10CM BAG/12

AWD13S



ADDITIONAL INFORMATION

Weight 0.2 kg

Dimensions $20 \times 25 \times 3 \text{ cm}$

Carton Qty 48

Pricing Unit Bag/12

Size 8cm x 10cm

Min Order Qty 1

ARTG ID 159433 **GMDN Code** 34970

Medical Device Class 1s

UNSPSC 42311545

Tax Schedule Taxable

Item Length 10

Item Width 8

Item Measurement Unit cm

Is Serialized Item No

Is Dangerous Good No

PRODUCT DESCRIPTION

The AeroWound Wound Dressings are a general purpose wound dressing for use in preventing infection and controlling bleeding. Available in sizes to suit wounds from finger injuries to large wounds.

Highly Absorbent

• Highly absorbent wound pad controls bleeding and absorbs wound exudates



AEROWOUND #13 WOUND DRESSING 8 X 10CM

AWD13S

Low adherent

nerent coating	minimizes w	ound trauma	associated wit	h changing	dressings
ł	herent coating	herent coating minimizes w	herent coating minimizes wound trauma	herent coating minimizes wound trauma associated wit	herent coating minimizes wound trauma associated with changing

Individually packed and sterile

• Dressings are available individually on demand in an easy tear wrapper and sterile for wound cleanliness.

Ease of application

• The wound pad is attached to a wrapping leader of gauze bandage making the dressing simple to apply, yet comfortable and secure for all general wound dressing.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Wound Dressing #13 Small