



AEROWOUND #13 WOUND DRESSING 8 X 10CM BAG/12

AWD13S

conforming bandage allows the pad to be positioned securely over contoured areas such as forearms, calves, and joints, helping the dressing stay in place during movement. Clear sizing across the range makes it straightforward to select the right coverage for fingertip injuries through to larger cuts.

Trusted for everyday first aid readiness, AeroWound Wound Dressing combines rapid application, reliable absorbency, and secure compression—ideal for kit builders and safety programs seeking consistent performance from a proven wound care staple.

- **Highly absorbent pad:** helps control bleeding and manages exudate.
- **Low-adherent contact layer:** supports gentle, low-pain removal.
- **Sterile:** promotes wound cleanliness at the time of application.
- **Integrated bandage leader:** enables firm, stable compression.
- **Hypoallergenic and latex-free:** designed for broad user comfort.

Sized **8 × 10 cm** for compact coverage with secure compression. Ideal for first response where a small, highly absorbent pad and integrated gauze leader provide quick placement and confident bleeding control.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Wound Dressing #13 Small