



## ADDITIONAL INFORMATION

<b>Weight</b>	0.25 kg
<b>Dimensions</b>	24 × 24 × 6.5 cm
<b>Carton Qty</b>	20
<b>Pricing Unit</b>	Bag/12
<b>Size</b>	18cm x 18cm
<b>Min Order Qty</b>	1
<b>ARTG ID</b>	159433
<b>GMDN Code</b>	34970
<b>Medical Device Class</b>	Class 1s
<b>UNSPSC</b>	42311545
<b>Tax Schedule</b>	Taxable
<b>Item Length</b>	18
<b>Item Width</b>	18
<b>Item Measurement Unit</b>	cm
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No
<b>Barcode (GTIN/EAN/UPC)</b>	9341394000711

## PRODUCT DESCRIPTION

The AeroWound Wound Dressings are a general purpose wound dressing for use in preventing infection and controlling bleeding. Available in sizes to suit wounds from finger injuries to large wounds.

### Highly Absorbent

- Highly absorbent wound pad controls bleeding and absorbs wound exudates



# AEROWOUND #15 WOUND DRESSING 18 X 18CM BAG/12

AWD15S

## **Low adherent**

- Special low adherent coating minimizes wound trauma associated with changing dressings

## **Individually packed and sterile**

- Dressings are available individually on demand in an easy tear wrapper and sterile for wound cleanliness.

## **Ease of application**

- The wound pad is attached to a wrapping leader of gauze bandage making the dressing simple to apply, yet comfortable and secure for all general wound dressing

Wound Dressing #15 Large