



## ADDITIONAL INFORMATION

<b>Weight</b>	0.25 kg
<b>Dimensions</b>	24 × 24 × 6.5 cm
<b>Carton Qty</b>	20
<b>Pricing Unit</b>	Bag/12
<b>Size</b>	18cm x 18cm
<b>Min Order Qty</b>	1
<b>ARTG ID</b>	159433
<b>GMDN Code</b>	34970
<b>Medical Device Class</b>	Class 1s
<b>UNSPSC Code</b>	4231 1545
<b>Tax Schedule</b>	Taxable
<b>Item Length</b>	18
<b>Item Width</b>	18
<b>Item Measurement Unit</b>	cm
<b>Is Medical Device</b>	Yes
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No



## PRODUCT DESCRIPTION

AeroWound Wound Dressing is a general-purpose compress designed to help control bleeding and protect a wide range of wounds. The soft, highly absorbent pad draws exudate away from the site while the low-adherent contact surface supports gentler removal. Built for dependable first aid, it provides rapid coverage and firm, stable compression when secured using the attached gauze bandage leader. Hypoallergenic materials support skin comfort for most users.

The dressing is sterile to support wound cleanliness at the point of care in workplaces, vehicles, schools, and home kits. The



# AEROWOUND #15 WOUND DRESSING 18 X 18CM BAG/12

AWD15S

conforming bandage allows the pad to be positioned securely over contoured areas such as forearms, calves, and joints, helping the dressing stay in place during movement. Clear sizing across the range makes it straightforward to select the right coverage for fingertip injuries through to larger cuts.

Trusted for everyday first aid readiness, AeroWound Wound Dressing combines rapid application, reliable absorbency, and secure compression—ideal for kit builders and safety programs seeking consistent performance from a proven wound care staple.

- **Highly absorbent pad:** helps control bleeding and manages exudate.
- **Low-adherent contact layer:** supports gentle, low-pain removal.
- **Sterile:** promotes wound cleanliness at the time of application.
- **Integrated bandage leader:** enables firm, stable compression.
- **Hypoallergenic and latex-free:** designed for broad user comfort.

Supplied as a single sterile compress measuring **18 cm × 18 cm** with an attached gauze bandage leader. The larger pad size offers generous coverage and secure compression for managing heavier bleeding during first response.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Wound Dressing #15 Large