



ADDITIONAL INFORMATION

Weight	0.0548 kg
Dimensions	15.2 × 2.6 × 2.6 cm
Carton Qty	100
Pricing Unit	Unit/1
Size	25g
Min Order Qty	10
ARTG ID	406566
GMDN Code	47764
Medical Device Class	Class 2b
Tax Schedule	Non-Taxable
Is Serialized Item	No
Is Dangerous Good	No
Barcode (GTIN/EAN/UPC)	9341394031159

PRODUCT DESCRIPTION

AeroWound Hydrogel offers healthcare professionals a trusted choice in advanced wound management. This amorphous wound care formulation has been designed for efficient adhesion and absorption. Fortified with propylene glycol, it adeptly balances moisture in the wound, attracting essential hydration while effectively managing excess exudates. AeroWound supports the body's natural autolytic debridement of devitalised tissue to promote faster wound healing.

<https://aerohealthcareonline.com/wp-content/uploads/2023/04/Amorphous-Gel-Final.mp4>

Features:

- Propylene Glycol Infusion:** Beyond just healing, AeroWound is all about creating an ideal environment for wounds. Enriched with propylene glycol, it actively draws in and seals moisture, even in drier wound settings.
- Optimal Moisture Management:** While ensuring hydration, AeroWound simultaneously absorbs excessive exudates, maintaining the perfect moisture balance for every wound type.



AEROWOUND AMORPHOUS WOUND HYDROGEL TUBE 25G

AWH25

3. **Autolytic Debridement Support:** AeroWound doesn't just treat the wound—it aids the body's natural processes. It facilitates the gentle autolytic debridement of devitalised tissue, clearing the way for healthier tissue regeneration.
4. **Accelerated Healing:** With the removal of non-viable tissue and optimal moisture levels, AeroWound sets the stage for quicker and more efficient wound healing.

Safety Considerations: For external use only. Always read the label and follow the directions for use.

Wound Hydrogel 25g