

AEROFIX NON-WOVEN ISLAND DRESSING 9 X 15CM BOX/3

AFID153S





ADDITIONAL INFORMATION

Weight 0.05 kg

Dimensions $1.5 \times 14 \times 20.3 \text{ cm}$

Carton Qty 96

Pricing Unit Box/3

Size 9 x 15cm

Min Order Qty 1

ARTG ID 159198 **GMDN Code** 44990

Medical Device Class 1s

UNSPSC 42311546

Tax Schedule Taxable

Item Length 15

Item Width 9

Item Measurement Unit cm

Is Medical Device Yes

Is Serialized Item No

Is Dangerous Good No

Shelf Life (Months) 36

PRODUCT DESCRIPTION

AEROFIX Non-Woven Island Dressing is a sterile, latex-free adhesive dressing designed for simple, effective coverage of larger wounds, abrasions and sutured sites. A soft, conformable non-woven backing helps the dressing sit smoothly over body contours while maintaining reliable fixation. The island design focuses coverage where it is needed and keeps edges low-profile for comfortable wear under clothing and during movement.

At the centre is the SupaSorb absorbent pad engineered to draw blood and exudate into internal pockets. This high-capacity pad supports a cleaner wound environment and reduces unnecessary changeovers. A non-adherent contact layer sits against the skin to minimise sticking, so removal is gentler and less disruptive to fragile or healing tissue.

Fixation is delivered by SupaBond adhesive, developed for strong, uniform adhesion that helps prevent lift and edge curl. The breathable, water-repellent non-woven barrier enhances patient comfort while providing dependable protection in everyday and



AEROFIX NON-WOVEN ISLAND DRESSING 9 X 15CM BOX/3

AFID153S

clinical settings. Individually sterile and intended for single use, AEROFIX offers clinicians, workplaces and home users a dependable choice for routine wound care.

- **Strong hold:** SupaBond adhesive helps prevent edge lift and displacement.
- **High absorbency:** SupaSorb pad manages exudate to keep the site cleaner.
- **Gentle removal:** non-adherent layer reduces trauma at dressing changes.
- Conformable comfort: soft non-woven backing contours to the body.
- **EO-sterile:** suitable for broad clinical and first aid use.

Size: 9 × 15 cm provides generous coverage for post-procedural sites and larger abrasions while remaining comfortable over curved areas, helping maintain a secure seal and reliable protection during movement.

Packaged in **boxes of 3** to support small procedures, top-ups and compact first aid kits, reducing waste while ensuring availability for larger wounds and sutured sites.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Aerofix non-woven Island Dressing