

AEROPLAST PREMIUM FABRIC FINGERTIP & KNUCKLE DRESSINGS BOX/12

AFP3012





ADDITIONAL INFORMATION

Weight 0.034 kg

Dimensions $3.2 \times 8.5 \times 10 \text{ cm}$

Carton Qty 120

Pricing Unit Box/12

Size Fingertip & Knuckle

Min Order Qty 6

ARTG ID 159198 **GMDN Code** 44990

Medical Device Class 1s

UNSPSC 42311546

Tax Schedule Taxable

Is Medical Device Yes
Is Serialized Item No
Is Dangerous Good No
Shelf Life (Months) 60

PRODUCT DESCRIPTION

The AeroPlast™ Premium Weight Fabric Adhesive Bandage is a top-tier, durable dressing designed for extended wear and maximum protection. Featuring our SupaBond™ extra-strong adhesive, this bandage is the ideal choice when long-lasting, reliable coverage is needed.

Aeroplast™ Premium Weight Fabric maximum long wearing durable, fabric for ultimate wound protection



AEROPLAST PREMIUM FABRIC FINGERTIP & KNUCKLE DRESSINGS BOX/12

AFP3012

Supabond™ Hypo-allergenic adhesive exclusively engineered, for low irritation and high adhesion

wound exudate pad super absorption, reducing skin maceration

Non-stick barrier protects wound with a non-adherent barrier, minimising trauma on removal

Unmatched SupaBond™ Adhesive

Engineered exclusively for AeroPlast™, the SupaBond™ adhesive offers unrivaled adhesion that secures firmly across the application site, preventing lift-up or edge curling. SupaBond™ is latex-free and leaves no sticky residue, addressing common issues found with heavy-duty fabric plasters.

Superior SupaSorb™ Wound Pad The SupaSorb™ wound pad is designed for optimal wound management, effectively absorbing blood and exudates into its internal pockets. This super-absorbent pad reduces the frequency of dressing changes, protecting the wound and keeping it clean and comfortable. The non-adherent film layer minimizes pain and trauma during dressing changes.

Heavy-Duty Fabric Backing

Constructed with conforming and air-permeable heavy-duty fabric, the AEROPLAST^M Premium Fabric dressing is perfect for demanding applications, such as in the manufacturing industry. Its superior adhesion and flexibility ensure the bandage stays in place and conforms to body contours, maintaining wound cleanliness and comfort under the toughest conditions. Knuckle and Fingertip Dressings are expertly crafted to provide superior protection and comfort for knuckle and fingertip injuries. Specifically shaped to fit the unique contours of knuckles and fingertips, providing optimal coverage and protection.

Available as a box of 12, this product is perfect for first aid kits, medical facilities, and industrial settings.

Contraindications and Precautions

Change dressing daily or more often according to the condition of the wound. As with all adhesive products apply and remove carefully from sensitive or fragile skin. Seek medical advice immediately for all serious wounds and burns. If redness or discomfort occurs, discontinue use and consult a physician.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

F/tip & Knuckle Dressings x 12



AEROPLAST PREMIUM FABRIC FINGERTIP & KNUCKLE DRESSINGS BOX/12

AFP3012