



ADDITIONAL INFORMATION

Weight	0.069 kg
Dimensions	6.5 × 6.7 × 10 cm
Carton Qty	48
Pricing Unit	Box/25
Size	7.5 x 4.5cm
Min Order Qty	12
ARTG ID	159198
GMDN Code	44990
Medical Device Class	Class 1s
UNSPSC	42311546
Tax Schedule	Taxable
Item Length	7.5
Item Width	4.5
Item Measurement Unit	cm
Is Serialized Item	No
Is Dangerous Good	No
Barcode (GTIN/EAN/UPC)	9341394005631

PRODUCT DESCRIPTION

The AeroPlast™ Premium Weight Fabric Adhesive Bandage is a top-tier, durable dressing designed for extended wear and maximum protection. Featuring our SupaBond™ extra-strong adhesive, this bandage is the ideal choice when long-lasting, reliable coverage is needed.



AEROPLAST PREMIUM FABRIC LARGE FINGERTIP 7.5 X 4.5CM BOX/25

AFP450

Aeroplast™
Premium Weight Fabric
maximum long wearing durable,
fabric for ultimate wound protection

Supabond™
Hypo-allergenic adhesive
exclusively engineered,
for low irritation and high adhesion

wound
exudate pad
super absorption,
reducing skin maceration

Non-stick barrier
protects wound with a non-adherent barrier,
minimising trauma on removal

Unmatched SupaBond™ Adhesive

Engineered exclusively for AeroPlast™, the SupaBond™ adhesive offers unrivaled adhesion that secures firmly across the application site, preventing lift-up or edge curling. SupaBond™ is latex-free and leaves no sticky residue, addressing common issues found with heavy-duty fabric plasters.

Superior SupaSorb™ Wound Pad The SupaSorb™ wound pad is designed for optimal wound management, effectively absorbing blood and exudates into its internal pockets. This super-absorbent pad reduces the frequency of dressing changes, protecting the wound and keeping it clean and comfortable. The non-adherent film layer minimizes pain and trauma during dressing changes.

Heavy-Duty Fabric Backing

Constructed with conforming and air-permeable heavy-duty fabric, the AEROPLAST™ Premium Fabric dressing is perfect for demanding applications, such as in the manufacturing industry. Its superior adhesion and flexibility ensure the bandage stays in place and conforms to body contours, maintaining wound cleanliness and comfort under the toughest conditions.

The fingertip 7.5cm x 4.5cm premium fabric bandage is designed for small to medium fingertip wounds. Offering a durable protective layer to create a moist healing environment.

Available as a box of 25

Contraindications and Precautions

Change dressing daily or more often according to the condition of the wound. As with all adhesive products apply and remove carefully from sensitive or fragile skin. Seek medical advice immediately for all serious wounds and burns. If redness or discomfort occurs, discontinue use and consult a physician.

This product is not a substitute for direct medical attention from a healthcare practitioner.

Fabric Standard Strip Box/100



**AEROPLAST PREMIUM FABRIC LARGE FINGERTIP
7.5 X 4.5CM BOX/25**

AFP450