



## ADDITIONAL INFORMATION

<b>Weight</b>	0.069 kg
<b>Dimensions</b>	6.5 × 6.7 × 10 cm
<b>Carton Qty</b>	48
<b>Pricing Unit</b>	Box/25
<b>Size</b>	7.5 x 4.5cm
<b>Min Order Qty</b>	12
<b>ARTG ID</b>	159198
<b>GMDN Code</b>	44990
<b>Medical Device Class</b>	Class 1s
<b>UNSPSC Code</b>	42311546
<b>Tax Schedule</b>	Taxable
<b>Item Length</b>	7.5
<b>Item Width</b>	4.5
<b>Item Measurement Unit</b>	cm
<b>Is Medical Device</b>	Yes
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No
<b>Shelf Life (Months)</b>	60
<b>Warranty Period (Months)</b>	12



## PRODUCT DESCRIPTION

AEROPLAST Premium Fabric adhesive bandages provide durable, long-wearing protection for minor cuts, grazes and abrasions in



# AEROPLAST PREMIUM FABRIC LARGE FINGERTIP 7.5 X 4.5CM BOX/25

AFP450

workplaces, education, hospitality and home first aid. A premium-weight woven backing conforms to fingers, knuckles and joints while allowing airflow to the skin for comfort. SupaBond extra-strong adhesive delivers high holding power that resists edge lift during frequent hand use and movement. The adhesive system is latex-free and hypo-allergenic for broad skin compatibility, supporting users with sensitive skin.

The dressing core features a SupaSorb absorbent pad with a non-adherent contact layer to manage light exudate and reduce trauma on removal. The breathable, open-weave fabric backing is air-permeable to help reduce moisture build-up under the dressing, while the conforming fit maintains coverage across contoured areas. The robust fabric resists abrasion in demanding environments and maintains adhesion under pressure, helping the plaster stay put for longer wear time.

Designed for reliability across busy settings, the premium construction combines strong adhesion, breathable comfort and skin-friendly materials. The range is suitable for first aid rooms, cabinets and mobile kits where consistent performance and fast application matter.

- **Strong adhesion:** SupaBond technology holds fast and resists edge curl.
- **Skin compatible:** hypo-allergenic, latex-free adhesive suitable for sensitive skin.
- **Breathable backing:** air-permeable fabric supports comfort during wear.
- **Absorbent, non-stick pad:** manages light exudate and eases removal.
- **Conforming fit:** flexible fabric maintains coverage over digits and joints.

**Size – Large Fingertip 7.5 × 4.5 cm:** contoured shape wraps securely around fingertips for focused coverage. Conforming, air-permeable fabric and strong adhesive help it stay in place. Skin-friendly and latex-free.

Supplied in a **box of 25** to streamline kit restocking and reduce downtime between replacements—ideal for workplace cabinets and service vehicles.

#### Contraindications and Precautions

Change dressing daily or more often according to the condition of the wound. As with all adhesive products apply and remove carefully from sensitive or fragile skin. Seek medical advice immediately for all serious wounds and burns. If redness or discomfort occurs, discontinue use and consult a physician.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Fabric Standard Strip Box/100