



ADDITIONAL INFORMATION

Weight	0.007 kg
Dimensions	11.5 × 7.6 × 1.9 cm
Carton Qty	25
Pricing Unit	Unit/1
Size	60mm
Min Order Qty	1
ARTG ID	257365
GMDN Code	42424
Medical Device Class	Class 1s
UNSPSC	42271913
Tax Schedule	Non-Taxable
Item Length	6
Item Measurement Unit	cm
Is Serialized Item	No
Is Dangerous Good	No

PRODUCT DESCRIPTION

Guedel Airways are a critical item for any advanced trauma response to keep airways clear of obstruction when providing advanced first response treatment.

Guedel Airways are available in a range of sizes which must be appropriately matched to the patient to ensure a correct fit. Colours help with ease of recognition of a suitable size. Each airway is individually packaged for best hygiene practice.

Guedel Airway #0 is colour black, and is 60mm in diameter.

Grey Guedel Airway Size 0



GREY CHILD STERILE GUEDEL AIRWAY SIZE 0 (ISO 5.5)

GAB1000