

GREEN SMALL STERILE GUEDEL AIRWAY SIZE 2 (ISO 8.0)

GAB1002



ADDITIONAL INFORMATION

Weight 0.013 kg

Dimensions $3.5 \times 11.6 \times 3.6 \text{ cm}$

Carton Qty 100
Pricing Unit Unit/1
Size 80mm

Min Order Qty 1

ARTG ID 257365
GMDN Code 42424
Medical Device Class 1s
UNSPSC 42271913

Tax Schedule Non-Taxable

Item Length 8

Item Measurement Unit cm

Is Medical Device Yes

Is Serialized Item No

Is Dangerous Good No

PRODUCT DESCRIPTION

Guedel Airways are a critical item for any advanced trauma response to keep airways clear of obstruction when providing advanced first response treatment.

Guedel Airways are available in a range of sizes which must be appropriately matched to the patient to ensure a correct fit. Colours help with ease of recognition of a suitable size. Each airway is individually packaged for best hygiene practice.

Guedel Airway #1 is colour green, and is 80mm in diameter.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.



GREEN SMALL STERILE GUEDEL AIRWAY SIZE 2 (ISO 8.0)

GAB1002

Green Guedel Airway Size 2