



## ADDITIONAL INFORMATION

<b>Weight</b>	0.02 kg
<b>Dimensions</b>	11.6 × 3.2 × 14.1 cm
<b>Carton Qty</b>	80
<b>Pricing Unit</b>	Unit/1
<b>Size</b>	90mm
<b>Min Order Qty</b>	1
<b>ARTG ID</b>	257365
<b>GMDN Code</b>	42424
<b>Medical Device Class</b>	Class 1s
<b>UNSPSC</b>	42271913
<b>Tax Schedule</b>	Non-Taxable
<b>Item Length</b>	9
<b>Item Measurement Unit</b>	cm
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No
<b>Barcode (GTIN/EAN/UPC)</b>	5030267079353

## PRODUCT DESCRIPTION

Guedel Airways are a critical item for any advanced trauma response to keep airways clear of obstruction when providing advanced first response treatment.

Guedel Airways are available in a range of sizes which must be appropriately matched to the patient to ensure a correct fit. Colours help with ease of recognition of a suitable size. Each airway is individually packaged for best hygiene practice.

Guedel Airway #3 is colour yellow, and is 90mm in diameter.



# YELLOW ADULT STERILE GUEDEL AIRWAY SIZE 3 (ISO 9.0)

GAB1003