

RED ADULT STERILE GUEDEL AIRWAY SIZE 4 (ISO

GAB1004



ADDITIONAL INFORMATION

Weight 0.02 kg

Dimensions $4.8 \times 11.6 \times 14.6 \text{ cm}$

Carton Qty 80

Pricing Unit Unit/1

Size 100mm

Min Order Qty 1

 ARTG ID
 257365

 GMDN Code
 42424

Medical Device Class 1s

UNSPSC 42271913

Tax Schedule Non-Taxable

Item Length10Item Measurement UnitcmIs Medical DeviceYesIs Serialized ItemNoIs Dangerous GoodNo

PRODUCT DESCRIPTION

Guedel Airways are a critical item for any advanced trauma response to keep airways clear of obstruction when providing advanced first response treatment.

Guedel Airways are available in a range of sizes which must be appropriately matched to the patient to ensure a correct fit. Colours help with ease of recognition of a suitable size. Each airway is individually packaged for best hygiene practice. Guedel Airway #4 is colour orange, and is 100mm in diameter.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.



RED ADULT STERILE GUEDEL AIRWAY SIZE 4 (ISO 10.0)

GAB1004