



ADDITIONAL INFORMATION

Weight	0.02 kg
Dimensions	4.8 × 11.6 × 14.6 cm
Carton Qty	80
Pricing Unit	Unit/1
Size	100mm
Min Order Qty	1
ARTG ID	257365
GMDN Code	42424
Medical Device Class	Class 1s
UNSPSC	42271913
Tax Schedule	Non-Taxable
Item Length	10
Item Measurement Unit	cm
Is Serialized Item	No
Is Dangerous Good	No
Barcode (GTIN/EAN/UPC)	5030267093533

PRODUCT DESCRIPTION

Guedel Airways are a critical item for any advanced trauma response to keep airways clear of obstruction when providing advanced first response treatment.

Guedel Airways are available in a range of sizes which must be appropriately matched to the patient to ensure a correct fit. Colours help with ease of recognition of a suitable size. Each airway is individually packaged for best hygiene practice.

Guedel Airway #4 is colour orange, and is 100mm in diameter.



RED ADULT STERILE GUEDEL AIRWAY SIZE 4 (ISO 10.0)

GAB1004