



ADDITIONAL INFORMATION

Weight	0.15 kg
Dimensions	8.5 × 8.5 × 10 cm
Carton Qty	50
Pricing Unit	Unit/1
Size	90cm x 11cm
Min Order Qty	1
ARTG ID	236793
GMDN Code	13681
Medical Device Class	Class 1
UNSPSC	42241507
Tax Schedule	Taxable
Item Length	90
Item Width	11
Item Measurement Unit	cm
Is Serialized Item	No
Is Dangerous Good	No
Barcode (GTIN/EAN/UPC)	19341394005003

PRODUCT DESCRIPTION

Light Weight Rolled Splint

Immobilisation is a critical item of care when responding to a casualty to limit the occurrence of further injury. Quality immobilisation using suitable products is integral to a full recovery of skeletal injuries.

Aluminium Foam Splints

An aluminium foam splint is a highly adaptive splint that can be rolled folded or formed into almost any application with the use of a quality support bandage. Available in a rolled or flat form, shrink-wrapped for protection.



AEROSPLINT ROLLED ALUMINIUM FOAM SPLINT 90 X 11CM

AR1018R

Aluminium Foam Splint Rolled