



## ADDITIONAL INFORMATION

<b>Weight</b>	0.15 kg
<b>Carton Qty</b>	50
<b>Pricing Unit</b>	Splint
<b>Size</b>	90cm x 11cm
<b>Min Order Qty</b>	1
<b>ARTG ID</b>	236793
<b>GMDN Code</b>	13681
<b>Medical Device Class</b>	Class 1
<b>UNSPSC</b>	42241507
<b>Tax Schedule</b>	Taxable
<b>Item Length</b>	90
<b>Item Width</b>	11
<b>Item Measurement Unit</b>	cm
<b>Location</b>	AU-Warehouse
<b>Is Medical Device</b>	Yes
<b>Is Serialized Item</b>	No

## PRODUCT DESCRIPTION

### Light Weight Flat Splint

Immobilisation is a critical item of care when responding to a casualty to limit the occurrence of further injury. Quality immobilisation using suitable products is integral to a full recovery of skeletal injuries.

### Aluminium Foam Splints

An aluminium foam splint is a highly adaptive splint that can be rolled folded or formed into almost any application with the use of a quality support bandage. Available in a rolled or flat form, shrink-wrapped for protection.

Aluminium Foam Splint Flat



# AEROSPLINT FOLDED ALUMINIUM FOAM 90 X 11CM

AR1018F