



## ADDITIONAL INFORMATION

<b>Weight</b>	0.04 kg
<b>Dimensions</b>	11 × 0.5 × 22.8 cm
<b>Carton Qty</b>	200
<b>Pricing Unit</b>	Unit/1
<b>Size</b>	22cm x 11cm
<b>Min Order Qty</b>	1
<b>ARTG ID</b>	236793
<b>GMDN Code</b>	13681
<b>Medical Device Class</b>	Class 1
<b>UNSPSC</b>	42241507
<b>Tax Schedule</b>	Taxable
<b>Item Length</b>	22
<b>Item Width</b>	11
<b>Item Measurement Unit</b>	cm
<b>Location</b>	AU-Warehouse
<b>Is Medical Device</b>	Yes
<b>Is Serialized Item</b>	No

## PRODUCT DESCRIPTION

### Light Weight Flat Splint

Immobilisation is a critical item of care when responding to a casualty to limit the occurrence of further injury. Quality immobilisation using suitable products is integral to a full recovery of skeletal injuries.

### Aluminium Foam Splints

An aluminium foam splint is a highly adaptive splint that can be rolled folded or formed into almost any application with the use of a quality support bandage. Available in a rolled or flat form, shrink-wrapped for protection.



# AEROSPLINT WRIST ALUMINIUM FOAM SPLINT 22 X 11CM

AR1019

Aluminium Foam Splint Flat