

AEROPLAST PLASTIC X-WIDE STRIP 7.5 X 2.5CM BOX/500

AP805





ADDITIONAL INFORMATION

Weight	0.42 kg
Dimensions	16 × 16.2 × 10.3 cm
Carton Qty	16
Pricing Unit	Box/500
Size	7.5 x 2.5cm
Min Order Qty	4
ARTG ID	159198
GMDN Code	44990
Medical Device Class	Class 1s
UNSPSC	42311546
Tax Schedule	Taxable
Item Length	7.5
ltem Width	2.5
Item Measurement Unit	cm
Is Medical Device	Yes
Is Serialized Item	No
Is Dangerous Good	No
Shelf Life (Months)	60

PRODUCT DESCRIPTION

The AEROPLAST Plastic range of adhesive plasters and bandages is designed to offer superior wound care solutions for everyday needs. These washproof and hypoallergenic bandages are crafted to provide comfort and protection for minor wounds, making them ideal for general applications. The innovative SupaSorb Wound Pad technology ensures excellent wound management by efficiently absorbing exudates, thereby reducing the frequency of dressing changes. This feature keeps the wound clean and comfortable, while the non-adherent film minimizes pain and trauma during dressing changes.

Equipped with the exclusive SupaBond adhesive, these bandages guarantee secure and lasting adhesion, preventing lift-up or edge curling, even in challenging conditions. The AEROPLAST Plastic bandages are made from materials that are both air permeable and conforming, allowing them to bend, stretch, and adapt to various shapes and forms. This ensures that the skin and wound can breathe as they heal, promoting a more comfortable recovery process.



AEROPLAST PLASTIC X-WIDE STRIP 7.5 X 2.5CM BOX/500

AP805

These bandages are latex-free, making them suitable for individuals with latex sensitivities. Manufactured to the highest standards, AEROPLAST Plastic bandages meet the needs of both medical professionals and everyday users. Whether you're a healthcare provider, a first aid trainer, or a family looking for reliable wound care solutions, the AEROPLAST Plastic range offers the protection and comfort you need. Trust in the quality and performance of these bandages to keep wounds protected and promote healing.

The AEROPLAST Plastic X-Wide Strip 7.5 x 2.5cm offers enhanced coverage for larger wounds, ensuring superior protection and comfort. Designed with the same washproof and hypoallergenic properties, this variant is ideal for general applications where extra width is required.

The AEROPLAST Plastic X-Wide Strip 7.5 x 2.5cm Box/500 provides a bulk supply of high-quality bandages, perfect for healthcare facilities and large households. This SKU ensures that you have ample bandages on hand for all your wound care needs.

Contraindications and Precautions

Change dressing daily or more often according to the condition of the wound. As with all adhesive products apply and remove carefully from sensitive or fragile skin. Seek medical advice immediately for all serious wounds and burns. If redness or discomfort occurs, discontinue use and consult a physician.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

AEROPLAST Extra Wide Strip