



## ADDITIONAL INFORMATION

<b>Pricing Unit</b>	pack/2
<b>Location</b>	AU-Warehouse
<b>Min Order Qty</b>	1
<b>ARTG ID</b>	227602
<b>GMDN Code</b>	34864
<b>Medical Device Class</b>	Class 1

## PRODUCT DESCRIPTION

The HyFin Vent Chest Seal Twin Pack is the industry leader in the treatment of penetrating injuries to the chest. Two separately packaged vented chest seals are offered in one unique package for treatment of both entry/exit or multiple penetrating injuries to the chest.

The HyFin® Vent Chest Seal design provides 3-vented channels that prevent airflow into the chest cavity during inhalation, while allowing air to escape through the vent channels during exhalation. The 3-vent channels allow blood to escape and also provide a backup fail-safe system, as even if two of the three channels become obstructed, the vent will remain fully operational.

Advanced adhesive technology provides superior adhesion in the most adverse conditions, including sweaty or hairy casualties. Packaged in rugged, easy-to-open foil pouches, the perforated packaging allows rescuers to open only one dressing at a time as needed. Each chest seal also includes a gauze pad to wipe the wound surface prior to application. Each HyFin® Vent Chest Seal has a large, pull tab for single-step peel-and-apply application and allows for the burping of the wound if necessary. The clear, transparent backing allows for easy placement over the wound area and conformability to the patient's chest.

HyFin Vent Chest Seal - Twin Pack