



ADDITIONAL INFORMATION

Weight	0.0266 kg
Dimensions	3.2 × 8.5 × 10 cm
Carton Qty	120
Pricing Unit	Box/50
Size	2.25 x 2.25cm
Min Order Qty	6
ARTG ID	159198
GMDN Code	44990
Medical Device Class	Class 1s
UNSPSC	42311546
Tax Schedule	Taxable
Item Length	2.25
Item Width	2.25
Item Measurement Unit	cm
Is Serialized Item	No
Is Dangerous Good	No
Barcode (GTIN/EAN/UPC)	9341394016330

PRODUCT DESCRIPTION

The AEROPLAST Transparent range brings to you a transparent aesthetic film blending with every skin tone and providing your wound with a waterproof barrier in a plaster form with a full range of sizes.

washproof sheer backing that bends and conforms to different shapes and forms

Supabond™
Hypo-allergenic adhesive
exclusively engineered,
for low irritation and high adhesion

SuperSorb™ Wound
Exudate Pad
super absorption,
reducing skin maceration

Non-stick Barrier
protects wound with a non-adherent barrier,
minimising trauma on removal

Unique transparent plaster backing

- The AeroPlast Transparent film is designed to bend, stretch and conform to body contours and guarantee freedom of movement. The film is aesthetic and waterproof, and its transparency assists in wound inspection.

SupaSorb™ Wound Pad

- An embossed wound pad provides superior wound management by draining blood or exudates into pockets in the wound pad. The Superabsorbent wound pad minimises the number of dressing changes whilst protecting the wound from further injury. SupaSorb™ leaves the wound site clean and comfortable, while the non-adherent film reduces the pain and site trauma when changing the dressing.

SupaBond™ Adhesive

- SupaBond adhesive, exclusively engineered for unrivalled adhesion, provides a secure adhesion across the site, preventing lift up or edge curling.

Comprehensive range of sizes



AEROPLAST TRANSPARENT SPOT DRESSINGS 2.25 X 2.25CM BOX/50

AT3035

- There is an AeroPlast Fabric dressing for any application with specialized fingertip and knuckle dressings, large and small patches, standards and extra wide plasters.

Contraindications / Precautions

Change dressing daily or more often according to the condition of the wound. As with all adhesive products apply and remove carefully from sensitive or fragile skin. Seek medical advice immediately for all serious wounds and burns. If redness or discomfort occurs, discontinue use and consult a physician.

This product is not a substitute for direct medical attention from a healthcare practitioner

Transparent Spot Dressing