

AEROPLAST TRANSPARENT SPOT DRESSINGS 2.25 X 2.25CM BOX/50

AT3035





ADDITIONAL INFORMATION

Weight 0.0266 kg

Dimensions $3.2 \times 8.5 \times 10 \text{ cm}$

Carton Qty 120

Pricing Unit Box/50

Size 2.25 x 2.25cm

Min Order Qty 6

ARTG ID 159198 **GMDN Code** 44990

Medical Device Class 1s

UNSPSC 42311546

Tax Schedule Taxable

Item Length 2.25

Item Width 2.25

Item Measurement Unit cm

Is Medical Device Yes

Is Serialized Item No

Is Dangerous Good No

Shelf Life (Months) 60

PRODUCT DESCRIPTION

The AEROPLAST Transparent range brings to you a transparent aesthetic film blending with every skin tone and providing your wound with a waterproof barrier in a plaster form with a full range of sizes.



AEROPLAST TRANSPARENT SPOT DRESSINGS 2.25 X 2.25CM BOX/50

AT3035

Aeroplast™ Plastic washproof sheer backing that bends and conforms to different shapes and forms

Supabond™ Hypo-allergenic adhesive exclusively engineered, for low irritation and high adhesion

SuperSorb™ Wound Exudate Pad super absorption, reducing skin maceration

Non-stick Barrier protects wound with a non-adherent barrier, minimising trauma on removal

Unique transparent plaster backing

• The AeroPlast Transparent film is designed to bend, stretch and conform to body contours and guarantee freedom of movement. The film is aesthetic and waterproof, and its transparency assists in wound inspection.

SupaSorb™ Wound Pad

An embossed wound pad provides superior wound management by draining blood or exudates into pockets in the wound pad.
 The Superabsorbent wound pad minimises the number of dressing changes whilst protecting the wound from further injury.
 SupaSorb™ leaves the wound site clean and comfortable, while the non-adherent film reduces the pain and site trauma when changing the dressing.

SupaBond™ Adhesive

• SupaBond adhesive, exclusively engineered for unrivalled adhesion, provides a secure adhesion across the site, preventing lift up or edge curling.

Comprehensive range of sizes

• There is an AeroPlast Fabric dressing for any application with specialized fingertip and knuckle dressings, large and small patches, standards and extra wide plasters.

Contraindications / Precautions

Change dressing daily or more often according to the condition of the wound. As with all adhesive products apply and remove carefully from sensitive or fragile skin. Seek medical advice immediately for all serious wounds and burns. If redness or discomfort occurs, discontinue use and consult a physician.

This product is not a substitute for direct medical attention from a healthcare practitioner

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.



AEROPLAST TRANSPARENT SPOT DRESSINGS 2.25 X 2.25CM BOX/50

AT3035

Transparent Spot Dressing