



ADDITIONAL INFORMATION

Weight	22.66 kg
Dimensions	130 × 23 × 67 cm
Carton Qty	1
Pricing Unit	Unit/1
Size	2200mm x 645mm x 80mm
Min Order Qty	1
UNSPSC	42171604
Tax Schedule	Taxable
Item Length	220
Item Width	64.5
Item Height	8
Item Measurement Unit	cm
Is Serialized Item	No
Is Dangerous Good	No
Barcode (GTIN/EAN/UPC)	29341394019564

PRODUCT DESCRIPTION

AeroRescue Stretchers are a light-weight high strength solution for many and varied patient recovery situations. Available in a variety of designs and materials for use in a wide range of environments and applications.

STR-09

Weight limit: 159kg

The collapsible basket stretcher has the same application as the STR-08 but detaches in half for ease of storage.

Folds to 126cm x 64.5cm x 21cm

Basket Stretcher Collapsible

