

ADDITIONAL INFORMATION

Weight	19 kg
Dimensions	225 × 65 × 20 cm
Carton Qty	1
Pricing Unit	Unit/1
Size	2200mm x 645mm x 165mm
Min Order Qty	1
UNSPSC	42171604
Tax Schedule	Taxable
Item Length	220
Item Width	64.5
Item Height	16.5
Item Measurement Unit	cm
Is Serialized Item	No
Is Dangerous Good	No
Barcode (GTIN/EAN/UPC)	29341394019557



PRODUCT DESCRIPTION

AeroRescue Stretchers are a light-weight high strength solution for many and varied patient recovery situations. Available in a variety of designs and materials for use in a wide range of environments and applications. The basket aviation stretcher is ideal for all types of aid in special situations such as lifting and transport by helicopter. Materials are fireproof and highly resistant to wear and corrosion.

Model: STR-08

Type: Basket

Applications: Emergency

Functioning: Lifting, Height-Adjustable

Other Characteristics: Plastic

Dimensions: 220 cm x 64.5 cm x 16.5 cm (86.6 in x 25.4 in x 6.5 in)

Weight: 18 kg (39.68 lb)

Loading Capacity: 159 kg (350.53 lb)

Additional Features:

- Suitable for a wide range of emergency situations
- Sturdy and flexible with quick reliable fittings and special sling equipment
- Ideal for all types of aid in special situations such as lifting and transport by helicopter
- Adjustable feet-securing mechanism, safety belts, and a mattress
- Designed to be secured inside the helicopter
- Materials are fireproof and highly resistant to wear and corrosion
- Does not release any toxic or polluting substances.

Basket Stretcher