



ADDITIONAL INFORMATION

Weight	0.05 kg
Dimensions	6.5 × 2 × 11 cm
Carton Qty	100
Pricing Unit	Unit/1
Size	65cm x 24cm
Min Order Qty	1
UNSPSC	42171703
Tax Schedule	Taxable
Item Length	65
Item Width	24
Item Measurement Unit	cm
Is Serialized Item	No
Is Dangerous Good	No
Barcode (GTIN/EAN/UPC)	9341394020825

PRODUCT DESCRIPTION

Hypothermia treatment is a critical part of first aid when dealing with trauma induced shock and where the body has insufficient ability to maintain body temperature.

Protect newborns from developing hypothermia immediately after birth with the Blizzard Baby Wrap™. This innovative 3-layer insulating pouch, made of Reflexcell material, delivers Tog 7 level of warmth that is both waterproof and windproof. A newborn may experience difficulty in regulating his/her core body temperature relative to the environmental temperature. The Blizzard Baby Wrap™ Reflexcell technology traps warm air, the metalized surface reflects heat back to the newborn, elasticated triple layer pouch draws material close to the body - reducing cold spaces & convection heat loss, thereby assisting in maintenance of a stable core body temperature. The Blizzard Baby Wrap™'s resealable double-sided tape closure easily allows paramedics to monitor the newborn while safeguarding the core body temperature. The Blizzard Baby Wrap™ is an important component of an EMS neonatal kit.

FEATURES

- 3-layer pouch
- Reflexcell material
- Double-sided tape closure
- Waterproof
- 7 TOG warmth
- Available color: Orange
- Flat-packed

Blizzard Baby Wrap