



ADDITIONAL INFORMATION

Weight	0.26 kg
Dimensions	29 × 19 × 1 cm
Carton Qty	25
Pricing Unit	Each
Size	233 x 196cm
Min Order Qty	1
Tax Schedule	Taxable
Item Length	233
Item Width	196
Item Measurement Unit	cm
Is Serialized Item	No
Is Dangerous Good	No
Barcode (GTIN/EAN/UPC)	5060297210259

PRODUCT DESCRIPTION

- Provides Tog 7 thermal insulation, preventing heat loss in the trauma patient in a range of cases, including water immersion, burns, inclement weather or hemorrhage
- Reflexcell® material's cellular network traps warm air, the metalized surface reflects heat to the patient, and the hem cinch cord prevents heat from escaping.
- Number 1 choice of the US Army
- The blanket is full-length with a 3-layer construction
- The Trauma Blanket is a reusable item and can be washed and kept ready for use in an emergency.

Blizzard Trauma Blanket



BLIZZARD GREEN TRAUMA BLANKET 233 X 196CM

BPS2-01