

BLUEY DISPOSABLE UNDERPAD 60 X 40CM

AMR335



ADDITIONAL INFORMATION

Weight 0.0244 kg

Dimensions $10.3 \times 8.4 \times 0.3 \text{ cm}$

Each

Carton Qty 300

Pricing Unit

Size 56 x 40cm

Min Order Qty 50

UNSPSC 42132101

Tax Schedule Non-Taxable

Item Length56Item Width40Item Measurement Unitcm

Is Medical Device No
Is Serialized Item No

Is Dangerous Good No

PRODUCT DESCRIPTION

Disposable Underpads are a highly absorbent and waterproof bed pad or mattress protector which is suitable for chair or bed protection. Also great as a disposable baby change mat or as an absorbent disposable pad for preparation of vaccines and other delicate medical procedures.

The Polypropylene backing layer protects the mattress and sheets by preventing fluid from passing through. The combination polymer and fluff core layer absorbs liquid away from the body and helps limit odour. Meanwhile the top non-woven layer is soft against the skin for maximum comfort.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.



BLUEY DISPOSABLE UNDERPAD 60 X 40CM

AMR335