



**ADDITIONAL INFORMATION**

<b>Weight</b>	0.0244 kg
<b>Dimensions</b>	10.3 × 8.4 × 0.3 cm
<b>Carton Qty</b>	250
<b>Pricing Unit</b>	Each
<b>Size</b>	56 x 40cm
<b>Min Order Qty</b>	1
<b>UNSPSC</b>	42132101
<b>Tax Schedule</b>	Non-Taxable
<b>Item Length</b>	56
<b>Item Width</b>	40
<b>Item Measurement Unit</b>	cm
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No

**PRODUCT DESCRIPTION**

Disposable Underpads are a highly absorbent and waterproof bed pad or mattress protector which is suitable for chair or bed protection. Also great as a disposable baby change mat or as an absorbent disposable pad for preparation of vaccines and other delicate medical procedures.

The Polypropylene backing layer protects the mattress and sheets by preventing fluid from passing through. The combination polymer and fluff core layer absorbs liquid away from the body and helps limit odour. Meanwhile the top non-woven layer is soft against the skin for maximum comfort.