



ADDITIONAL INFORMATION

Weight	4.5 kg
Dimensions	62 × 19.5 × 32.4 cm
Carton Qty	1
Pricing Unit	Each
Size	Baby
Min Order Qty	1
UNSPSC Code	42301502
Tax Schedule	Taxable
Is Serialized Item	No
Is Dangerous Good	No



PRODUCT DESCRIPTION

The Brayden Baby is the new arrival to 'our family' range of Brayden manikins. Brayden Baby has been designed to be an ideal teaching manikin for Paediatric BLS, EPALS and EPILS. Developed with help from leading European experts in the field of infant resuscitation, Brayden Baby CPR Manikin's unique features include intuitive lights that provide real-time feedback for ventilation training and lights providing real-time feedback for chest compressions (depth, rate, recoil and finger/thumb position). When all component parts of chest compressions are performed correctly, a CPR quality indication light illuminates on the forehead of the Brayden Baby manikin. Both sets of lights help guide the student to perform good quality infant CPR.

Aero Healthcare is the exclusive Australian importer of the entire [Brayden range](#), including the [Brayden Baby](#), the [Brayden Pro](#), [Brayden Pulse](#), [Brayden Junior](#) and the [Brayden](#).

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

