



## ADDITIONAL INFORMATION

<b>Weight</b>	4.64 kg
<b>Dimensions</b>	61 × 33.6 × 25.9 cm
<b>Carton Qty</b>	1
<b>Pricing Unit</b>	Unit/1
<b>Size</b>	Adult
<b>Min Order Qty</b>	1
<b>UNSPSC</b>	42301502
<b>Tax Schedule</b>	Taxable
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No



## PRODUCT DESCRIPTION

The Brayden CPR Manikin lets your students see the real point of CPR getting blood flowing to the brain and other organs with LED lights that illuminate the blood flow from the heart to the brain during CPR. The lights can be switched on or off as required.

Winner of the 2016 iF Design Award for Medicine & Health, this essential visual learning tool instantly reinforces CPR Trainers instructions to students. The Brayden CPR Manikin 's dual modes help structure CPR training: from learning effective CPR to perfecting CPR technique, in accordance with ARC guidelines.

[Click Here for Our Dedicated Brayden Site.](#)

Package Contents Includes:

- 10 x Face Masks
- 1 x Extra Lung

- 6 x Batteries
- 1 x Face Skin
- 1 x Instructions
- 1 x Manikin Carry Bag

## Technical Specifications

Model:	IM13-R	IM16	IM13-S
Type:	Torso	Torso	Torso
Display visual flow of blood from the heart to the brain:	Yes	Yes	No
Bluetooth Connectivity / Sensor Data Collection:	No	Yes	No
Chest compression depth and speed measuring sensor:	Yes	Yes	No
Chest movement with ventilation:	Yes	Yes (when unpaired with App)	Yes
Head tilt and Chin lift for opening airway:	Yes	Yes	Yes
Airways system:	Non-breathing with disposable lungs	Non-breathing with disposable lungs	Non-breathing with disposable lungs
Ribs, Xiphoid:	Yes	Yes	Yes
Audible feedback:	Clicker on/off	Clicker on/off	Clicker on/off
Length:	59cm	59cm	59cm
Weight:	6.8lb (including batteries)	7lb (including batteries)	6.6lb
Operating temperature:	59°F to 86°F	32°F to 104°F	59°F to 86°F
Storage temperature:	14°F to 104°F	-4°F to 140°F	14°F to 104°F
Humidity:	Relative humidity between 5% and 95%	Relative humidity between 5% and 95%	Relative humidity between 5% and 95%
Warranty:	2 years (excluding consumables)	2 years (excluding consumables)	2 years (excluding consumables)
Power:	AC Adapter (Input AC 100 ~ 260V, Output DC 9V 1A) Battery C type (6 batteries)	AC Adapter (Input AC 100 ~ 260V, Output DC 9V 1A) Battery C type (6 batteries)	No

Brayden CPR Manikin (Red Blood Circulation)