



## ADDITIONAL INFORMATION

|                           |                     |
|---------------------------|---------------------|
| <b>Weight</b>             | 4.64 kg             |
| <b>Dimensions</b>         | 61 × 33.6 × 25.9 cm |
| <b>Carton Qty</b>         | 1                   |
| <b>Pricing Unit</b>       | Unit/1              |
| <b>Size</b>               | Adult               |
| <b>Min Order Qty</b>      | 1                   |
| <b>UNSPSC Code</b>        | 42301502            |
| <b>Tax Schedule</b>       | Taxable             |
| <b>Is Serialized Item</b> | No                  |
| <b>Is Dangerous Good</b>  | No                  |



## PRODUCT DESCRIPTION

The Brayden CPR Manikin lets your students see the real point of CPR getting blood flowing to the brain and other organs with LED lights that illuminate the blood flow from the heart to the brain during CPR. The lights can be switched on or off as required.

Winner of the 2016 iF Design Award for Medicine & Health, this essential visual learning tool instantly reinforces CPR Trainers instructions to students. The Brayden CPR Manikin 's dual modes help structure CPR training: from learning effective CPR to perfecting CPR technique, in accordance with ARC guidelines.

[Click Here for Our Dedicated Brayden Site.](#)

Package Contents Includes:

- 10 x Face Masks

- 1 x Extra Lung
- 6 x Batteries
- 1 x Face Skin
- 1 x Instructions
- 1 x Manikin Carry Bag

## Technical Specifications

| Model:   | IM13-R  | IM16  | IM13-S                               |
|--|---|---|--------------------------------------|
| <b>Type:</b>   | Torso   | Torso   | Torso                                |
| <b>Display visual flow of blood from the heart to the brain:</b> | Yes   | Yes   | No                                   |
| <b>Bluetooth Connectivity / Sensor Data Collection:</b>          | No  | Yes   | No                                   |
| <b>Chest compression depth and speed measuring sensor:</b>       | Yes   | Yes   | No                                   |
| <b>Chest movement with ventilation:</b>                          | Yes   | Yes (when unpaired with App)  | Yes                                  |
| <b>Head tilt and Chin lift for opening airway:</b>               | Yes   | Yes   | Yes                                  |
| <b>Airways system:</b>   | Non-breathing with disposable lungs   | Non-breathing with disposable lungs   | Non-breathing with disposable lungs  |
| <b>Ribs, Xiphoid:</b>  | Yes   | Yes   | Yes                                  |
| <b>Audible feedback:</b>   | Clicker on/off  | Clicker on/off  | Clicker on/off                       |
| <b>Length:</b>   | 59cm  | 59cm  | 59cm                                 |
| <b>Weight:</b>   | 6.8lb (including batteries)   | 7lb (including batteries)   | 6.6lb                                |
| <b>Operating temperature:</b>                                    | 59°F to 86°F  | 32°F to 104°F   | 59°F to 86°F                         |
| <b>Storage temperature:</b>                                      | 14°F to 104°F   | -4°F to 140°F   | 14°F to 104°F                        |
| <b>Humidity:</b>   | Relative humidity between 5% and 95%  | Relative humidity between 5% and 95%  | Relative humidity between 5% and 95% |
| <b>Warranty:</b>   | 2 years (excluding consumables)   | 2 years (excluding consumables)   | 2 years (excluding consumables)      |
| <b>Power:</b>  | AC Adapter (Input AC 100 ~ 260V, Output DC 9V 1A)<br>Battery C type (6 batteries) | AC Adapter (Input AC 100 ~ 260V, Output DC 9V 1A)<br>Battery C type (6 batteries) | No                                   |

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Brayden CPR Manikin (Red Blood Circulation)