

BRAYDEN PULSE CPR MANIKIN (WITH LIGHTS)

IM13-R



ADDITIONAL INFORMATION

Weight 4.64 kg

Dimensions $61 \times 33.6 \times 25.9 \text{ cm}$

Carton Qty 1

Pricing Unit Unit/1
Size Adult

Min Order Qty 1

UNSPSC 42301502
Tax Schedule Taxable

Is Serialized Item No
Is Dangerous Good No

PRODUCT DESCRIPTION

The Brayden CPR Manikin lets your students see the real point of CPR getting blood flowing to the brain and other organs with LED lights that illuminate the blood flow from the heart to the brain during CPR. The lights can be switched on or off as required.

Winner of the 2016 iF Design Award for Medicine & Health, this essential visual learning tool instantly reinforces CPR Trainers instructions to students. The Brayden CPR Manikin 's dual modes help structure CPR training: from learning effective CPR to perfecting CPR technique, in accordance with ARC guidelines.

Click Here for Our Dedicated Brayden Site.

Package Contents Includes:

- 10 x Face Masks
- 1 x Extra Lung



BRAYDEN PULSE CPR MANIKIN (WITH LIGHTS)

IM13-R

- 6 x Batteries
- 1 x Face Skin
- 1 x Instructions
- 1 x Manikin Carry Bag

Technical Specifications

Model:	IM13-R	IM16	IM13-S
Туре:	Torso	Torso	Torso
Display visual flow of blood			
from the heart to the brain:	Yes	Yes	No
Bluetooth Connectivity /			
Sensor Data Collection:	No	Yes	No
Chest compression depth			
·	Yes	Yes	No
and speed measuring sensor:	163	103	140
Chest movement with			
ventilation:	Yes	Yes (when unpaired with App)	Yes
Head tilt and Chin lift for	Yes	Yes	Yes
opening airway:			
Airways system:	Non-breathing with disposable lungs	Non-breathing with disposable lungs	Non-breathing with disposable lungs
Ribs, Xiphoid:	Yes	Yes	Yes
Audible feedback:	Clicker on/off	Clicker on/off	Clicker on/offs
Length:	59cm	59cm	59cm
Weight:	6.8lb (including batteries)	7lb (including batteries)	6.6lb
Operating temperature:	59°F to 86°F	32°F to 104°F	59°F to 86°F
Storage temperature:	14°F to 104°F	-4°F to 140°F	14°F to 104°F
Humidity:	Relative humidity between 5% and 95%	Relative humidity between 5% and 95%	Relative humidity between 5% and 95%
Warranty:	2 years (excluding consumables)	2 years (excluding consumables)	2 years (excluding consumables)
Power:	AC Adapter (Input AC 100 ~ 260V, Output DC 9V 1A)	AC Adapter (Input AC 100 ~ 260V, Output DC 9V 1A)	No
	Battery C type (6 batteries)	Battery C type (6 batteries)	

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.



BRAYDEN PULSE CPR MANIKIN (WITH LIGHTS)

IM13-R

Brayden CPR Manikin (Red Blood Circulation)