





### **ADDITIONAL INFORMATION**

Weight 4.06 kg

**Dimensions**  $61 \times 33.6 \times 25.9 \text{ cm}$ 

Carton Qty 1

Pricing Unit Unit/1

**Size** Adult

Min Order Qty 1

**UNSPSC** 42301502

**Tax Schedule** Taxable

**Item Length** 59

Item Width 33

**Item Height** 20

**Item Measurement Unit** cm

Is Serialized Item No

Is Dangerous Good No

### **PRODUCT DESCRIPTION**

With all the ventilation, airflow, realistic skin and audible clicker features of the IM13-R Brayden CPR Manikin, this manikin has no LED light function.

Package Contents Includes:

- 10 x Face Masks
- 1 x Extra Lung
- 1 x Face Skin
- 1 x Instructions
- 1 x Manikin Carry Bag



## **BRAYDEN CPR MANIKIN (WITHOUT LIGHTS)**

IM13-S

Also available in the Brayden <u>Pulse</u> and <u>Pro</u>.

For our dedicated Brayden Manikin Site, Click HERE.

### **Technical Specifications**

Model:	IM13-R	IM16	IM13-S
Туре:	Torso	Torso	Torso
Display visual flow of blood			
	Yes	Yes	No
from the heart to the brain:			
Bluetooth Connectivity / Sensor Data Collection:	No	Yes	No
Chest compression depth			
chest compression depth			
and speed measuring	Yes	Yes	No
sensor:			
Chest movement with	Yes	Yes (when unpaired with App)	Yes
ventilation:	163	res (when dispance with App)	163
Head tilt and Chin lift for	Yes	Yes	Yes
opening airway:	No. 1	No. 1 detail the St. Personal I.	No. of the state o
Airways system:	Non-breathing with disposable lungs	Non-breathing with disposable lungs	Non-breathing with disposable lungs
Ribs, Xiphoid:	Yes	Yes	Yes
Audible feedback:	Clicker on/off	Clicker on/off	Clicker on/offs
Length:	59cm	59cm	59cm
Weight:	6.8lb (including batteries)	7lb (including batteries)	6.6lb
Operating temperature:	59°F to 86°F	32°F to 104°F	59°F to 86°F
Storage temperature:	14°F to 104°F	-4°F to 140°F	14°F to 104°F
Humidity:	Relative humidity between 5%	Relative humidity between 5%	Relative humidity between
	and 95%	and 95%	5% and 95%
Warranty:	2 years (excluding	2 years (excluding consumables)	2 years (excluding
-	consumables)		consumables)
	AC Adapter (Input AC 100 ~ 260V, Output DC 9V 1A)	AC Adapter (Input AC 100 ~ 260V, Output DC 9V 1A)	
Power:	2004, Output DC 34 1A)	200V, Output DC 3V TA)	No
	Battery C type (6 batteries)	Battery C type (6 batteries)	



# **BRAYDEN CPR MANIKIN (WITHOUT LIGHTS)**

IM13-S