




ADDITIONAL INFORMATION

Weight	4.06 kg
Dimensions	61 × 33.6 × 25.9 cm
Carton Qty	1
Pricing Unit	Unit/1
Size	Adult
Min Order Qty	1
UNSPSC	42301502
Tax Schedule	Taxable
Item Length	59
Item Width	33
Item Height	20
Item Measurement Unit	cm
Is Serialized Item	No
Is Dangerous Good	No
Barcode (GTIN/EAN/UPC)	9341394015371



PRODUCT DESCRIPTION

With all the ventilation, airflow, realistic skin and audible clicker features of the IM13-R Brayden CPR Manikin, this manikin has no LED light function.

Package Contents Includes:

- 10 x Face Masks
- 1 x Extra Lung
- 1 x Face Skin
- 1 x Instructions

- 1 x Manikin Carry Bag

Also available in the Brayden [Pulse](#) and [Pro](#).

For our dedicated Brayden Manikin Site, [Click HERE](#).

Technical Specifications

Model:	IM13-R	IM16	IM13-S
Type:	Torso	Torso	Torso
Display visual flow of blood from the heart to the brain:	Yes	Yes	No
Bluetooth Connectivity / Sensor Data Collection:	No	Yes	No
Chest compression depth and speed measuring sensor:	Yes	Yes	No
Chest movement with ventilation:	Yes	Yes (when unpaired with App)	Yes
Head tilt and Chin lift for opening airway:	Yes	Yes	Yes
Airways system:	Non-breathing with disposable lungs	Non-breathing with disposable lungs	Non-breathing with disposable lungs
Ribs, Xiphoid:	Yes	Yes	Yes
Audible feedback:	Clicker on/off	Clicker on/off	Clicker on/off
Length:	59cm	59cm	59cm
Weight:	6.8lb (including batteries)	7lb (including batteries)	6.6lb
Operating temperature:	59°F to 86°F	32°F to 104°F	59°F to 86°F
Storage temperature:	14°F to 104°F	-4°F to 140°F	14°F to 104°F
Humidity:	Relative humidity between 5% and 95%	Relative humidity between 5% and 95%	Relative humidity between 5% and 95%
Warranty:	2 years (excluding consumables)	2 years (excluding consumables)	2 years (excluding consumables)
Power:	AC Adapter (Input AC 100 ~ 260V, Output DC 9V 1A) Battery C type (6 batteries)	AC Adapter (Input AC 100 ~ 260V, Output DC 9V 1A) Battery C type (6 batteries)	No

