



ADDITIONAL INFORMATION

Weight	4.06 kg
Dimensions	61 × 33.6 × 25.9 cm
Carton Qty	1
Pricing Unit	Unit/1
Size	Adult
Min Order Qty	1
UNSPSC Code	42301502
Tax Schedule	Taxable
Item Length	59
Item Width	33
Item Height	20
Item Measurement Unit	cm
Is Serialized Item	No
Is Dangerous Good	No



PRODUCT DESCRIPTION

With all the ventilation, airflow, realistic skin and audible clicker features of the IM13-R Brayden CPR Manikin, this manikin has no LED light function.

Package Contents Includes:

- 10 x Face Masks
- 1 x Extra Lung

- 1 x Face Skin
- 1 x Instructions
- 1 x Manikin Carry Bag

Also available in the Brayden [Pulse](#) and [Pro](#).

For our dedicated Brayden Manikin Site, [Click HERE](#).

Technical Specifications

Model:	IM13-R	IM16	IM13-S
Type:	Torso	Torso	Torso
Display visual flow of blood from the heart to the brain:	Yes	Yes	No
Bluetooth Connectivity / Sensor Data Collection:	No	Yes	No
Chest compression depth and speed measuring sensor:	Yes	Yes	No
Chest movement with ventilation:	Yes	Yes (when unpaired with App)	Yes
Head tilt and Chin lift for opening airway:	Yes	Yes	Yes
Airways system:	Non-breathing with disposable lungs	Non-breathing with disposable lungs	Non-breathing with disposable lungs
Ribs, Xiphoid:	Yes	Yes	Yes
Audible feedback:	Clicker on/off	Clicker on/off	Clicker on/off
Length:	59cm	59cm	59cm
Weight:	6.8lb (including batteries)	7lb (including batteries)	6.6lb
Operating temperature:	59°F to 86°F	32°F to 104°F	59°F to 86°F
Storage temperature:	14°F to 104°F	-4°F to 140°F	14°F to 104°F
Humidity:	Relative humidity between 5% and 95%	Relative humidity between 5% and 95%	Relative humidity between 5% and 95%
Warranty:	2 years (excluding consumables)	2 years (excluding consumables)	2 years (excluding consumables)
Power:	AC Adapter (Input AC 100 ~ 260V, Output DC 9V 1A) Battery C type (6 batteries)	AC Adapter (Input AC 100 ~ 260V, Output DC 9V 1A) Battery C type (6 batteries)	No

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Brayden CPR Manikin (no illuminators)