



**ADDITIONAL INFORMATION**

<b>Weight</b>	4.5 kg
<b>Carton Qty</b>	1
<b>Pricing Unit</b>	Each
<b>Min Order Qty</b>	1
<b>Tax Schedule</b>	Taxable
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No



**PRODUCT DESCRIPTION**

Brayden Baby Pro resembles the Brayden Baby Advanced in appearance and shares many consumable components. Provides unique intuitive real-time ventilation feedback using lights and sound. The Brayden Baby Pro also facilitates comprehensive objective/numeric information of infant CPR performance when connected via Bluetooth with the [Brayden Online iOS app for iPad](#). Its flexibility in use with the Brayden Online product range and added features for scenario use makes it the ideal choice for training infant CPR in both hospital and pre-hospital settings.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

