



ADDITIONAL INFORMATION

Weight	0.92 kg
Dimensions	43 × 32 × 15 cm
Pricing Unit	Each
Min Order Qty	1
Tax Schedule	Taxable
Is Serialized Item	No
Is Dangerous Good	No

PRODUCT DESCRIPTION

Help reduce healthcare inequalities one heartbeat at a time

- Ethnically Accurate - Nasal and mouth features of an Ethnically African adult male to give a realistic ventilation experience
- Cost effective - A single manikin with multiple skins can provide diverse training experiences
- Light Weight & Portable - Easy and quick conversion from **Brayden Adult** to **Brayden OBI**
- Feedback lights - Brayden OBI skin work just as other Brayden skins

Tone/Texture: Medium Taupe/Matte

Easy AED training pad removal

Latex-free and durable material

Includes:

- 1 x Adult OBI body skin.
- 1 x Adult OBI face skin.

To be used on the Brayden Adult Basic/Advanced/Pro manikins.

BRAYDEN OBI Conversion Kit (compatible with IM13, IM16)