

BRAYDEN PRO CPR MANIKIN (BLUETOOTH CONNECTIVITY)

IM16-R





ADDITIONAL INFORMATION

Weight	4.62 kg	
Dimensions	61 × 33.6 × 25.9 cm	
Carton Qty	1	
Pricing Unit	Unit/1	
Size	Adult	
Min Order Qty	1	
UNSPSC	42301502	
Tax Schedule	Taxable	
Is Serialized Item	No	
Is Dangerous Good	No	

PRODUCT DESCRIPTION

All the benefits of the IM13-R Brayden CPR Manikin, plus sensor data output via bluetooth to the free Brayden Pro App. This provides trainers quantitative data both in real time and for review to improve the quality of CPR performance and provide a numeric assessment.

The Brayden Pro provides healthcare professionals with quantitative data-driven feedback.

The Brayden Pro sends live data to the <u>Brayden Pro App</u> via Bluetooth to provide real-time feedback of CPR performance this data can then be reviewed and analysed to improve the quality of CPR performance and provide a numeric, quantitative assessment. Up to 6 Brayden Pro 's can be connected to 1 device allowing for simultaneous training and assessment. Records can be stored and easily exported via Email for easy analysis on your Spreadsheet program.

Package Contents Includes:

• 10 x Face Masks



- 1 x Extra Lung
- 6 x Batteries
- 1 x Face Skin
- 1 x Instructions
- 1 x Manikin Carry Bag

Click Here for our Dedicated Brayden Site.

Technical Specifications

Model:	IM13-R	IM16	IM13-S
Туре:	Torso	Torso	Torso
Display visual flow of blood			
from the heart to the brain:	Yes	Yes	No
Bluetooth Connectivity /			
Sensor Data Collection:	No	Yes	No
Chest compression depth			
and speed measuring	Yes	Yes	No
sensor:			
Chest movement with ventilation:	Yes	Yes (when unpaired with App)	Yes
Head tilt and Chin lift for opening airway:	Yes	Yes	Yes
Airways system:	Non-breathing with disposable lungs	Non-breathing with disposable lungs	Non-breathing with disposable lungs
Ribs, Xiphoid:	Yes	Yes	Yes
Audible feedback:	Clicker on/off	Clicker on/off	Clicker on/offs
Length:	59cm	59cm	59cm
Weight:	6.8lb (including batteries)	7lb (including batteries)	6.6lb
Operating temperature:	59°F to 86°F	32°F to 104°F	59°F to 86°F
Storage temperature:	14°F to 104°F	-4°F to 140°F	14°F to 104°F
Humidity:	Relative humidity between 5% and 95%	Relative humidity between 5% and 95%	Relative humidity between 5% and 95%
Warranty:	2 years (excluding consumables)	2 years (excluding consumables)	2 years (excluding consumables)

IM16-R



BRAYDEN PRO CPR MANIKIN (BLUETOOTH CONNECTIVITY)

IM16-R

Power:

AC Adapter (Input AC 100 ~ 260V, Output DC 9V 1A)

AC Adapter (Input AC 100 ~ 260V, Output DC 9V 1A)

No

Battery C type (6 batteries)

Battery C type (6 batteries)

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Brayden PRO Manikin (Red Blood Circulation)