



## ADDITIONAL INFORMATION

<b>Weight</b>	4.62 kg
<b>Dimensions</b>	61 × 33.6 × 25.9 cm
<b>Carton Qty</b>	1
<b>Pricing Unit</b>	Unit/1
<b>Size</b>	Adult
<b>Min Order Qty</b>	1
<b>UNSPSC Code</b>	42301502
<b>Tax Schedule</b>	Taxable
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No



## PRODUCT DESCRIPTION

All the benefits of the IM13-R Brayden CPR Manikin, plus sensor data output via bluetooth to the free Brayden Pro App. This provides trainers quantitative data both in real time and for review to improve the quality of CPR performance and provide a numeric assessment.

The **Brayden Pro** provides healthcare professionals with **quantitative data-driven feedback**.

The Brayden Pro sends live data to the [Brayden Pro App](#) via Bluetooth to provide real-time feedback of CPR performance this data can then be reviewed and analysed to improve the quality of CPR performance and provide a numeric, quantitative assessment. Up to 6 Brayden Pro 's can be connected to 1 device allowing for simultaneous training and assessment. Records can be stored and easily exported via Email for easy analysis on your Spreadsheet program.

Package Contents Includes:

- 10 x Face Masks
- 1 x Extra Lung

- 6 x Batteries
- 1 x Face Skin
- 1 x Instructions
- 1 x Manikin Carry Bag

[Click Here for our Dedicated Brayden Site.](#)

## Technical Specifications

Model:	IM13-R	IM16	IM13-S
<b>Type:</b>	Torso	Torso	Torso
<b>Display visual flow of blood from the heart to the brain:</b>	Yes	Yes	No
<b>Bluetooth Connectivity / Sensor Data Collection:</b>	No	Yes	No
<b>Chest compression depth and speed measuring sensor:</b>	Yes	Yes	No
<b>Chest movement with ventilation:</b>	Yes	Yes (when unpaired with App)	Yes
<b>Head tilt and Chin lift for opening airway:</b>	Yes	Yes	Yes
<b>Airways system:</b>	Non-breathing with disposable lungs	Non-breathing with disposable lungs	Non-breathing with disposable lungs
<b>Ribs, Xiphoid:</b>	Yes	Yes	Yes
<b>Audible feedback:</b>	Clicker on/off	Clicker on/off	Clicker on/off
<b>Length:</b>	59cm	59cm	59cm
<b>Weight:</b>	6.8lb (including batteries)	7lb (including batteries)	6.6lb
<b>Operating temperature:</b>	59°F to 86°F	32°F to 104°F	59°F to 86°F
<b>Storage temperature:</b>	14°F to 104°F	-4°F to 140°F	14°F to 104°F
<b>Humidity:</b>	Relative humidity between 5% and 95%	Relative humidity between 5% and 95%	Relative humidity between 5% and 95%
<b>Warranty:</b>	2 years (excluding consumables)	2 years (excluding consumables)	2 years (excluding consumables)
<b>Power:</b>	AC Adapter (Input AC 100 ~ 260V, Output DC 9V 1A) Battery C type (6 batteries)	AC Adapter (Input AC 100 ~ 260V, Output DC 9V 1A) Battery C type (6 batteries)	No

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the



# BRAYDEN PRO CPR MANIKIN (BLUETOOTH CONNECTIVITY)

IM16-R

directions for use.

Brayden PRO Manikin (Red Blood Circulation)