



ADDITIONAL INFORMATION

Weight	0.14 kg
Dimensions	22 × 2 × 27.5 cm
Carton Qty	1
Pricing Unit	Each
Min Order Qty	1
Tax Schedule	Non-Taxable
Is Medical Device	Yes
Is Serialized Item	No
Is Dangerous Good	No

PRODUCT DESCRIPTION

The Cardiac Science Powerheart G3 Adult Defibrillation Pads (non-polarised) are designed for use with all Powerheart Automated External Defibrillators excluding the Powerheart G5 AED. The pads are non-polarised, meaning they can be placed in either position on a patient's torso, reducing confusion and pressure in a stressful scenario.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.