



## ADDITIONAL INFORMATION

<b>Weight</b>	0.46 kg
<b>Dimensions</b>	16 × 7 × 4 cm
<b>Carton Qty</b>	1
<b>Pricing Unit</b>	Each
<b>Min Order Qty</b>	1
<b>Tax Schedule</b>	Non-Taxable
<b>Is Serialized Item</b>	Yes
<b>Is Dangerous Good</b>	Yes

## PRODUCT DESCRIPTION

The G5 AED is a single, medical-grade Li-ion battery for use with the Cardiac Science Powerheart G5 AED. The battery can deliver up to 500 shocks at 150VE. 4 year expiry.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.