



ADDITIONAL INFORMATION

Weight	0.14 kg
Dimensions	2 × 27.5 cm
Carton Qty	1
Pricing Unit	Each
Min Order Qty	1
Tax Schedule	Non-Taxable
Is Serialized Item	No
Is Dangerous Good	No

PRODUCT DESCRIPTION

The G5 AED is a single, medical-grade Li-ion battery for use with the Cardiac Science Powerheart G5 AED. The battery can deliver up to 500 shocks at 150VE. 4 year expiry

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.